



食品分析

3-3

34 种 C6~C24 饱和/不饱和脂肪酸甲酯的分析

组份: C6:0、C8:0、C10:0、C11:0、C12:0、C13:0、C14:0、C14:1 n-5、C15:0、C15:1 n-5、C16:0、C16:1 n-7
 C17:0、C17:1 n-7、C18:0、C18:1 cis n-9、C18:1 trans n-9、C18:2 n-6、C18:3 n-6、C18:3 n-3、
 C18:4 n-3、C20:0、C20:1 n-9、C20:2 n-6、C20:3 n-6、C20:4 n-6、C20:3 n-3、C20:5、C22:0、
 C22:1 n-9、C22:2 n-6、C24:0、C22:6 n-3、C24:1 n-9 脂肪酸甲酯

色谱柱: MEGA-WAX 25m x 0.32mm 0.25um df

订货号: S-WAX-032-025-25

